

PACO49778

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Transcription repressor capable of inhibiting the transactivation capability of TCF3/E47. May play a role in regulating antigen-dependent B-cell differentiation.

Gene ID:

MSC

Uniprot

O60682

Synonyms:

Musculin (Activated B-cell factor 1) (ABF-1) (Class A basic helix-loop-helix protein 22) (bHLHa22), MSC, ABF1 BHLHA22


Immunogen:

Recombinant Human Musculin protein (3-143AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human lymph node tissue using PACO49778 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO49778 at dilution of 1:100.