

PROP1 Antibody, FITC conjugated



PACO49916

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Possibly involved in the ontogenesis of pituitary gonadotropes, as well as somatotropes, lactotropes and caudomedial thyrotropes.

Gene ID:

PROP1

Uniprot

O75360

Synonyms:

Homeobox protein prophet of Pit-1 (PROP-1) (Pituitary-specific homeodomain factor), PROP1

Immunogen:

Recombinant Human Homeobox protein prophet of Pit-1 protein (6-214AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A