

RPL12 Antibody



PACO49966

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Binds directly to 26S ribosomal RNA.

Gene ID:

RPL12

Uniprot

P30050

Synonyms:

60S ribosomal protein L12 (Large ribosomal subunit protein uL11), RPL12


Immunogen:

Recombinant Human 60S ribosomal protein L12 protein (1-99AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human liver tissue using PACO49966 at dilution of 1:100.