

RPS2 Antibody, FITC conjugated



PACO49988

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cell-cell adherens junction, cytoplasm, cytosol, cytosolic small ribosomal subunit, extracellular exosome, focal adhesion, membrane, nucleoplasm, nucleus, cadherin binding involved in cell-cell adhesion.

Gene ID:

RPS2

Uniprot

P15880

Synonyms:

40S ribosomal protein S2 (40S ribosomal protein S4) (Protein LLRep3) (Small ribosomal subunit protein uS5), RPS2, RPS4

Immunogen:

Recombinant Human 40S ribosomal protein S2 protein (2-61AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A