

# SKIL Antibody, FITC conjugated



PACO50004

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May have regulatory role in cell division or differentiation in response to extracellular signals.

**Gene ID:**

SKIL

**Uniprot**

P12757

**Synonyms:**

Ski-like protein (Ski-related oncogene) (Ski-related protein), SKIL, SNO

**Immunogen:**

Recombinant Human Ski-like protein (359-579AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A