

# SOX4 Antibody



PACO50026

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:1000-1:5000,  
IHC:1:20-1:200

**Protein Background:**

Transcriptional activator that binds with high affinity to the T-cell enhancer motif 5'-AACAAAG-3' motif.

**Gene ID:**

SOX4

**Uniprot**

Q06945

**Synonyms:**

Transcription factor SOX-4, SOX4

**Immunogen:**


Recombinant Human Transcription factor SOX-4 protein (1-60AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human brain tissue using PACO50026 at dilution of 1:100.



Western Blot. Positive WB detected in: HeLa whole cell lysate, 293 whole cell lysate, HepG2 whole cell lysate. All lanes: SOX4 antibody at 1:2000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 48 kDa. Observed band size: 48 kDa.