

STAT5B Antibody, FITC conjugated



PACO50044

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Carries out a dual function: signal transduction and activation of transcription. Mediates cellular responses to the cytokine KITLG/SCF and other growth factors. Binds to the GAS element and activates PRL-induced transcription.

Gene ID:

STAT5B

Uniprot

P51692

Synonyms:

Signal transducer and activator of transcription 5B, STAT5B

Immunogen:

Recombinant Human Signal transducer and activator of transcription 5B protein (1-787AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A