

# STAT5B Antibody, Biotin conjugated



PACO50045

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Carries out a dual function: signal transduction and activation of transcription. Mediates cellular responses to the cytokine KITLG/SCF and other growth factors. Binds to the GAS element and activates PRL-induced transcription.

**Gene ID:**

STAT5B

**Uniprot**

P51692

**Synonyms:**

Signal transducer and activator of transcription 5B, STAT5B

**Immunogen:**

Recombinant Human Signal transducer and activator of transcription 5B protein (1-787AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A