

TPST1 Antibody, Biotin conjugated



PACO50077

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the O-sulfation of tyrosine residues within acid, c motifs of polypeptides.

Gene ID:

TPST1

Uniprot

O60507

Synonyms:

Protein-tyrosine sulfotransferase 1 (EC 2.8.2.20) (Tyrosylprotein sulfotransferase 1) (TPST-1), TPST1

Immunogen:

Recombinant Human Protein-tyrosine sulfotransferase 1 protein (283-370AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A