

DUOXA1 Antibody



PACO50114

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

May be required for the maturation and the transport from the endoplasmic reticulum to the plasma membrane of functional DUOX1.

Gene ID:

DUOXA1

Uniprot

Q1HG43

Synonyms:

Dual oxidase maturation factor 1 (Dual oxidase activator 1) (Numb-interacting protein), DUOXA1, NIP NUMBIP


Immunogen:

Recombinant Human Dual oxidase maturation factor 1 protein (73-183AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunofluorescent analysis of HeLa cells using PACO50114 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).