

# HPX Antibody, FITC conjugated



PACO50176

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Binds heme and transports it to the liver for breakdown and iron recovery, after which the free hemopexin returns to the circulation.

**Gene ID:**

HPX

**Uniprot**

P02790

**Synonyms:**

Hemopexin (Beta-1B-glycoprotein), HPX

**Immunogen:**

Recombinant Human Hemopexin protein (32-254AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A