

cyp1a Antibody



PACO50222

Product Information

Size:

50ug

Reactivity:

Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Cytochromes P450 are a group of heme-thiolate monooxygenases.

Gene ID:

cyp1a

Uniprot

Q8UW07

Synonyms:

Immunogen:

Recombinant Zebrafish Cytochrome P450 protein (1-519AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A