

nkx2.5 Antibody



PACO50226

Product Information

Size:

50ug

Reactivity:

Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cardiac chamber development, cardiac muscle cell differentiation, cardiac muscle cell proliferation, cardiac ventricle development, embryonic heart tube development, heart development, heart looping, heart morphogenesis, pharyngeal arch artery morphogenesis, regulation of cardiac chamber formation.

Gene ID:

nkx2.5

Uniprot

Q567L4

Synonyms:

Immunogen:

Recombinant Zebrafish Nkx2.5 protein (1-312AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A