

CDH-1 Antibody, HRP conjugated



PACO50255

Product Information

Size:

50ug

Reactivity:

Phanerochaete chrysosporium

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Degrades both lignin and cellulose. Oxidizes cellobiose to cellobionolactone.

Gene ID:

CDH-1

Uniprot

Q01738

Synonyms:

Cellobiose dehydrogenase (CDH) (EC 1.1.99.18) (Cellobiose-quinone oxidoreductase), CDH-1; CDH-2

Immunogen:

Recombinant Phanerochaete chrysosporium Cellobiose dehydrogenase protein (19-208AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A