

Defensin-like protein 1 Antibody, Biotin conjugated

PACO50265



Product Information

Size:

50ug

Reactivity:

Dahlia merckii

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Possesses antimicrobial activity sensitive to inorganic cations. Has no inhibitory effect on insect gut alpha-amylase. Induces potential changes in fungal membranes and increased K⁺ efflux and Ca²⁺ uptake. Interacts with sphingolipids and ergosterols found in fungal plasma membranes.

Gene ID:**Uniprot**

P0C8Y4

Synonyms:

Defensin-like protein 1 (Cysteine-rich antimicrobial protein 1) (Defensin AMP1) (DmAMP1)

Immunogen:

Recombinant Dahlia merckii Defensin-like protein 1 protein (1-50AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A