

Sporamin A Antibody



PACO50278

Product Information

Size:

50ug

Reactivity:

Ipomoea batatas

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Major tuberous root protein.

Gene ID:

Uniprot

P14715

Synonyms:

Sporamin A (Clone PIM0335)

Immunogen:

Recombinant Ipomoea batatas Sporamin A protein (24-219AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A