

D8L Antibody



PACO50290

Product Information

Size:

50ug

Reactivity:

Vaccinia virus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Binds to chondroitin sulfate on the cell surface to provide virion attachment to target cell.

Gene ID:

D8L

Uniprot

P20508

Synonyms:

Cell surface-binding protein (Carbonic anhydrase homolog)


Immunogen:

Recombinant Vaccinia virus Cell surface-binding protein (1-304AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: D8L antibody at 3 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 52 kDa. Observed band size: 52 kDa.