

PACO50306

Product Information

Size:

50ug

Reactivity:

Bacillus subtilis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Probably involved in peptidoglycan modification during cortex synthesis.

Gene ID:

pbpE

Uniprot

P32959

Synonyms:

Penicillin-binding protein 4* (PBP 4*) (PBP 4A) (Penicillin-binding protein E), pbpE


Immunogen:

Recombinant Bacillus subtilis Penicillin-binding protein 4* protein (1-451AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: pbpE antibody at 3 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 52 kDa. Observed band size: 52 kDa.