

PPX1 Antibody, Biotin conjugated



PACO50313

Product Information

Size:

50ug

Reactivity:

Saccharomyces cerevisiae

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Degradation of inorganic polyphosphates.

Gene ID:

PPX1

Uniprot

P38698

Synonyms:

Exopolyphosphatase (ExopolyPase) (EC 3.6.1.11) (Metaphosphatase), PPX1

Immunogen:

Recombinant Saccharomyces cerevisiae Exopolyphosphatase protein (1-397AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A