

AFP2 Antibody



PACO50326

Product Information

Size:

50ug

Reactivity:

Raphanus sativus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Possesses antifungal activity sensitive to inorganic cations. Induces potential changes in fungal membranes and increased K(+) efflux and Ca(2+) uptake.

Gene ID:

AFP2

Uniprot

P30230

Synonyms:

Defensin-like protein 2 (Cysteine-rich antifungal protein 2) (AFP2) (RAFP2), AFP2


Immunogen:

Recombinant Raphanus sativus Defensin-like protein 2 protein (30-80AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: AFP2 antibody at 3 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. predicted band size: 22 kDa. observed band size: 22 kDa.