

pepF Antibody



PACO50330

Product Information

Size:

50ug

Reactivity:

Mycoplasma pneumoniae

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Gene ID:

pepF

Uniprot

P54125

Synonyms:

Oligoendopeptidase F homolog (EC 3.4.24), pepF


Immunogen:

Recombinant Mycoplasma pneumoniae Oligoendopeptidase F homolog protein (1-210AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: pepF antibody at 2.8 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. predicted band size: 40 kDa. observed band size: 40 kDa.