

ompK Antibody



PACO50338

Product Information

Size:

50ug

Reactivity:

Vibrio parahaemolyticus serotype O3: K6

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Serves as receptor for a broad-host-range vibriophage, KVP40.

Gene ID:

ompK

Uniprot

P59570

Synonyms:

Outer membrane protein OmpK, ompK


Immunogen:

Recombinant Vibrio parahaemolyticus serotype O3: K6 Outer membrane protein OmpK protein (21-266AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: ompK antibody at 3.4 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. predicted band size: 44 kDa. observed band size: 44 kDa.