

# hlgB Antibody, FITC conjugated



PACO50376

---

## Product Information

**Size:**

50ug

**Reactivity:**

Staphylococcus aureus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Toxin that seems to act by forming pores in the membrane of the cell. Has a hemolytic and a leucotoxic activity.

**Gene ID:**

hlgB

**Uniprot**

P0A075

**Synonyms:**

Gamma-hemolysin component B (H-gamma-1) (H-gamma-I), hlgB

**Immunogen:**

Recombinant Staphylococcus aureus Gamma-hemolysin component B protein (26-325AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A