

hlgB Antibody, Biotin conjugated



PACO50377

Product Information

Size:

50ug

Reactivity:

Staphylococcus aureus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Toxin that seems to act by forming pores in the membrane of the cell. Has a hemolytic and a leucotoxic activity.

Gene ID:

hlgB

Uniprot

P0A075

Synonyms:

Gamma-hemolysin component B (H-gamma-1) (H-gamma-I), hlgB

Immunogen:

Recombinant Staphylococcus aureus Gamma-hemolysin component B protein (26-325AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A