

RNU2 Antibody



PACO50382

Product Information

Size:

50ug

Reactivity:

Ustilago sphaerogena

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Gene ID:

RNU2

Uniprot

P00654

Synonyms:

Ribonuclease U2 (RNase U2) (EC 3.1.27.4), RNU2


Immunogen:

Recombinant Ustilago sphaerogena Ribonuclease U2 protein (1-114AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: RNU2 antibody at 3.25 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. predicted band size: 30 kDa. observed band size: 30 kDa.