

# gpd-1 Antibody



PACO50386

---

## Product Information

**Size:**

50ug

**Reactivity:**

Caenorhabditis elegans

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:****Gene ID:**

gpd-1

**Uniprot**

P04970

**Synonyms:**

Glyceraldehyde-3-phosphate dehydrogenase 1 (GAPDH-1) (EC 1.2.1.12), gpd-1

**Immunogen:**

Recombinant Caenorhabditis elegans Glyceraldehyde-3-phosphate dehydrogenase 1 protein (1-341AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A