

gpd-1 Antibody, HRP conjugated



PACO50387

Product Information

Size:

50ug

Reactivity:

Caenorhabditis elegans

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

gpd-1

Uniprot

P04970

Synonyms:

Glyceraldehyde-3-phosphate dehydrogenase 1 (GAPDH-1) (EC 1.2.1.12), gpd-1

Immunogen:

Recombinant Caenorhabditis elegans Glyceraldehyde-3-phosphate dehydrogenase 1 protein (1-341AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A