

# GOLGA3 Antibody, Biotin conjugated



PACO50461

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Golgi auto-antigen; probably involved in maintaining Golgi structure.

**Gene ID:**

GOLGA3

**Uniprot**

Q08378

**Synonyms:**

Golgin subfamily A member 3 (Golgi complex-associated protein of 170 kDa) (GCP170) (Golgin-160), GOLGA3

**Immunogen:**

Recombinant Human Golgin subfamily A member 3 protein (217-405AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A