

# PTK6 Antibody, HRP conjugated



PACO50499

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PTK6

**Uniprot**

Q13882

**Synonyms:**

Protein-tyrosine kinase 6 (EC 2.7.10.2) (Breast tumor kinase) (Tyrosine-protein kinase BRK), PTK6, BRK

**Immunogen:**

Recombinant Human Protein-tyrosine kinase 6 protein (1-113AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A