

PTK6 Antibody, Biotin conjugated



PACO50501

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

PTK6

Uniprot

Q13882

Synonyms:

Protein-tyrosine kinase 6 (EC 2.7.10.2) (Breast tumor kinase) (Tyrosine-protein kinase BRK), PTK6, BRK

Immunogen:

Recombinant Human Protein-tyrosine kinase 6 protein (1-113AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A