

SMN1 Antibody, FITC conjugated



PACO50512

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

SMN1

Uniprot

Q16637

Synonyms:

Survival motor neuron protein (Component of gems 1) (Gemin-1), SMN1; SMN2, SMN SMNT; SMNC

Immunogen:

Recombinant Human Survival motor neuron protein (3-194AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A