

FOXF2 Antibody



PACO50590

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probable transcription activator for a number of lung-specific genes.

Gene ID:

FOXF2

Uniprot

Q12947

Synonyms:

Forkhead box protein F2 (Forkhead-related activator 2) (FREAC-2) (Forkhead-related protein FKHL6) (Forkhead-related transcription factor 2), FOXF2, FKHL6 FREAC2

Immunogen:

Recombinant Human Forkhead box protein F2 protein (226-364AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A