

PACO50606

---

## Product Information

**Size:**

50ug

**Reactivity:**

Brucella abortus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Transcriptional regulator mediating the effects of the quorum sensing autoinducer C12-HSL (N-dodecanoyl-homoserine lactone).

**Gene ID:**

vjbR

**Uniprot**

Q2YJ50

**Synonyms:**

HTH-type quorum sensing-dependent transcriptional regulator VjbR, vjbR

**Immunogen:**

Recombinant Brucella abortus HTH-type quorum sensing-dependent transcriptional regulator VjbR protein (1-259AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A