

**Guanine nucleotide-binding protein subunit
beta-like protein Antibody, HRP conjugated
PACO50611**



Product Information

Size:

50ug

Reactivity:

Glycine max

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Uniprot

Q39836

Synonyms:

Guanine nucleotide-binding protein subunit beta-like protein

Immunogen:

Recombinant Glycine max Guanine nucleotide-binding protein subunit beta-like protein (1-325AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A