

vacA Antibody



PACO50614

Product Information

Size:

50ug

Reactivity:

Helicobacter pylori

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Induces vacuolation of eukaryotic cells. Causes ulceration and gastric lesions.

Gene ID:

vacA

Uniprot

Q48253

Synonyms:

Vacuolating cytotoxin autotransporter [Cleaved into: Vacuolating cytotoxin; Vacuolating cytotoxin translocator], vacA

Immunogen:

Recombinant Helicobacter pylori Vacuolating cytotoxin autotransporter protein (1038-1310AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A