

FHL5 Antibody, FITC conjugated



PACO50652

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in the regulation of spermatogenesis. Stimulates CREM transcriptional activity in a phosphorylation-independent manner.

Gene ID:

FHL5

Uniprot

Q5TD97

Synonyms:

Four and a half LIM domains protein 5 (FHL-5) (Activator of cAMP-responsive element modulator in testis) (Activator of CREM in testis), FHL5, ACT

Immunogen:

Recombinant Human Four and a half LIM domains protein 5 protein (16-284AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A