

SZT2 Antibody



PACO50654

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Involved in oxidative stress. May be involved in superoxide dismutase activity and in neuroprotective effect of peroxisomes, but has no direct dismutase activity when tested in yeast. May enhance epileptogenesis and confer low seizure threshold.

Gene ID:

SZT2

Uniprot

Q5T011

Synonyms:

KICSTOR complex protein SZT2 (Seizure threshold 2 protein homolog), SZT2, C1orf84 KIAA0467


Immunogen:

Recombinant Human Protein SZT2 protein (2396-2517AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO50654 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO50654 at dilution of 1:100.