

TCTN3 Antibody, Biotin conjugated



PACO50669

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Part of the tectonic-like complex which is required for tissue-specific ciliogenesis and may regulate ciliary membrane composition. May be involved in apoptosis regulation. Necessary for signal transduction through the sonic hedgehog (Shh) signaling pathway.

Gene ID:

TCTN3

Uniprot

Q6NUS6

Synonyms:

Tectonic-3, TCTN3, C10orf61 TECT3

Immunogen:

Recombinant Human Tectonic-3 protein (348-472AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A