

# CRTC1 Antibody, Biotin conjugated



PACO50697

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CRTC1

**Uniprot**

Q6UUU9

**Synonyms:**

CREB-regulated transcription coactivator 1 (Mucoepidermoid carcinoma translocated protein 1) (Transducer of regulated cAMP response element-binding protein 1) (TORC-1) (Transducer of CREB protein 1), CRTC1, KIAA0616 MECT1 TORC1 WAMTP1

**Immunogen:**

Recombinant Human CREB-regulated transcription coactivator 1 protein (1-248AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A