

# TRIM59 Antibody



PACO50754

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

May serve as a multifunctional regulator for innate immune signaling pathways.

**Gene ID:**

TRIM59

**Uniprot**

Q8IWR1

**Synonyms:**

Tripartite motif-containing protein 59 (RING finger protein 104) (Tumor suppressor TSBF-1), TRIM59, RNF104 TRIM57 TSBF1

**Immunogen:**


Recombinant Human Tripartite motif-containing protein 59 protein (126-181AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human liver cancer using PACO50754 at dilution of 1:100.