

TRIM59 Antibody, Biotin conjugated



PACO50757

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May serve as a multifunctional regulator for innate immune signaling pathways.

Gene ID:

TRIM59

Uniprot

Q8IWR1

Synonyms:

Tripartite motif-containing protein 59 (RING finger protein 104) (Tumor suppressor TSBF-1), TRIM59, RNF104 TRIM57 TSBF1

Immunogen:

Recombinant Human Tripartite motif-containing protein 59 protein (126-181AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A