

# UNC45B Antibody



PACO50774

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Acts as a co-chaperone for HSP90 and is required for proper folding of the myosin motor domain. Plays a role in sarcomere formation during muscle cell development. Is necessary for normal early lens development.

**Gene ID:**

UNC45B

**Uniprot**

Q8IWX7

**Synonyms:**

Protein unc-45 homolog B (Unc-45B) (SMUNC45), UNC45B, CMYA4 UNC45

**Immunogen:**


Recombinant Human Protein unc-45 homolog B protein (602-710AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO50774 at dilution of 1:100.