

LOXHD1 Antibody



PACO50778

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Involved in hearing. Required for normal function of hair cells in the inner ear.

Gene ID:

LOXHD1

Uniprot

Q8IVV2

Synonyms:

Lipoxygenase homology domain-containing protein 1, LOXHD1


Immunogen:

Recombinant Human Lipoxygenase homology domain-containing protein 1 protein (1101-1366AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human pancreatic cancer using PACO50778 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO50778 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).