

# DYNC2H1 Antibody, Biotin conjugated



PACO50785

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May function as a motor for intraflagellar retrograde transport. Functions in cilia biogenesis. May play a role in transport between endoplasmic reticulum and Golgi or organization of the Golgi in cells.

**Gene ID:**

DYNC2H1

**Uniprot**

Q8NCM8

**Synonyms:**

Cytoplasmic dynein 2 heavy chain 1 (Cytoplasmic dynein 2 heavy chain) (Dynein cytoplasmic heavy chain 2) (Dynein heavy chain 11) (hDHC11) (Dynein heavy chain isotype 1B), DYNC2H1, DHC1B DHC2 DNCH2 DYH1B KIAA1997

**Immunogen:**

Recombinant Human Cytoplasmic dynein 2 heavy chain 1 protein (1928-2065AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A