

# PF0142 Antibody



PACO50810

---

## Product Information

**Size:**

50ug

**Reactivity:**

Pyrococcus furiosus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:5000

**Protein Background:**

**Gene ID:**

PF0142

**Uniprot**

Q8U4E6

**Synonyms:**

Putative L-asparaginase (EC 3.5.1.1) (L-asparagine amidohydrolase) [Cleaved into: Putative L-asparaginase subunit alpha; Putative L-asparaginase subunit beta]

**Immunogen:**


Recombinant Pyrococcus furiosus Putative L-asparaginase protein (1-175AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Western Blot. Positive WB detected in Recombinant protein. All lanes: PF0142 antibody at 2.4 $\mu$ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 36 kDa. Observed band size: 36 kDa.