

# APOBEC4 Antibody, Biotin conjugated



PACO50817

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Putative C to U editing enzyme whose physiological substrate is not yet known.

**Gene ID:**

APOBEC4

**Uniprot**

Q8WW27

**Synonyms:**

Putative C->U-editing enzyme APOBEC-4 (EC 3.5.4) (Apolipoprotein B mRNA-editing enzyme catalytic polypeptide-like 4), APOBEC4, C1orf169

**Immunogen:**

Recombinant Human Putative C->U-editing enzyme APOBEC-4 protein (221-367AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A