

DRAM1 Antibody



PACO50830

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Lysosomal modulator of autophagy that plays a central role in p53/TP53-mediated apoptosis. Not involved in p73/TP73-mediated autophagy.

Gene ID:

DRAM1

Uniprot

Q8N682

Synonyms:

DNA damage-regulated autophagy modulator protein 1 (Damage-regulated autophagy modulator), DRAM1, DRAM


Immunogen:

Recombinant Human DNA damage-regulated autophagy modulator protein 1 protein (50-100AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human lung cancer using PACO50830 at dilution of 1:100.