

COQ8B Antibody, Biotin conjugated



PACO50841

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in CoQ10 (COQ10A and/or COQ10B) biosynthesis, which is required for podocyte migration.

Gene ID:

COQ8B

Uniprot

Q96D53

Synonyms:

Atypical kinase COQ8B, mitochondrial (EC 2.7) (AarF domain-containing protein kinase 4) (Coenzyme Q protein 8B), COQ8B, ADCK4

Immunogen:

Recombinant Human Atypical kinase COQ8B, mitochondrial protein (10-137AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A