

# BBX Antibody, FITC conjugated



PACO50868

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Transcription factor that is necessary for cell cycle progression from G1 to S phase.

**Gene ID:**

BBX

**Uniprot**

Q8WY36

**Synonyms:**

HMG box transcription factor BBX (Bobby sox homolog) (HMG box-containing protein 2), BBX, HBP2

**Immunogen:**

Recombinant Human HMG box transcription factor BBX protein (62-292AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A