

FAM111A Antibody



PACO50914

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Chromatin-associated protein required for PCNA loading on replication sites. Promotes S-phase entry and DNA synthesis. May directly function at replication forks, explaining why Simian virus 40 (SV40) interacts with FAM111A to overcome host range restriction.

Gene ID:

FAM111A

Uniprot

Q96PZ2

Synonyms:

Protein FAM111A, FAM111A, KIAA1895


Immunogen:

Recombinant Human Protein FAM111A protein (26-179AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human liver cancer using PACO50914 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO50914 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).